

Rules for bright field and DIC microscopy:

1. Köhler illumination
  - Switch on the light source
  - Focus sample and then move out from sample area
  - Close down field stop to get octagonal shape
  - Move condenser to see sharp edges of this octagonal shape
  - Put this octagonal shape in a middle of view area
  - Open field stop far enough to not see black colour
2. Set the light
  - Try to use 3200K for pictures
  - Decrease the light intensity by 2 filter rings sitting on right side of microscope (close to micro/macro screw)
  - Check aperture setting
  - Choose technique you want to use on condenser ring (bright field - H/I, DIC – I, II or III)
  - Check settings on filter ring over objective, in case of DIC choose position 8, in case of bright field choose position 7
3. Epifluorescence:
  - Switch on source of fluorescent light
  - Choose correct filter according to list below
  - **Do not forget to switch light off!!!**
4. Troubleshooting
  - If you see yellow / green / blue light in the microscope, check filter ring over objectives
  - If view is disturb by some dark object, check all filters position, most probably it is not left “in a middle position”
  - If you see your sample correctly in a microscope but on a computer it looks yellow/ pinkish, play with white balance in camera settings
  - If you see some stable, old “burn out” sample or any view on a monitor which you don’t see in a microscope – contact Marta

| Microscope position | Filter set |   | Usage  |
|---------------------|------------|---|--|
| <b>1</b>            | #DAPI      | Exc. BP 365<br>B. splitter FT 395<br>Em. LP 397           | Fluorescence<br><b>DAPI</b>                        |
| <b>2</b>            | #9         | Exc. BP 450-490<br>B. splitter FT 510<br>Em. LP 515       | Fluorescence<br><b>GFP/YFP</b><br>broadband filter |
| <b>3</b>            | #15        | Exc. BP 546+/-12<br>B. splitter FT 580<br>Em. LP 590      | Fluorescence<br><b>RFP</b><br>broadband filter     |
| <b>4</b>            | #46        | Exc. BP 500+/-20<br>B. splitter FT 515<br>Em. LP 535+/-30 | Fluorescence<br><b>GFP, Alexa 488, FITC etc.</b>   |
| <b>5</b>            | #47        | Exc. BP 436+/-20<br>B. splitter FT 455<br>Em. LP 480+/-40 | Fluorescence<br><b>CFP</b>                         |
| <b>6</b>            | #38        | Exc. BP 470+/-40<br>B. splitter FT 515<br>Em. LP 525+/-50 | Fluorescence<br><b>CFP+GFP</b>                     |
| <b>7</b>            | empty      |   | <b>Bright field pic.</b>                           |
| <b>8</b>            | DIC        | Analyzer (polarization filter DIC)                        | <b>DIC</b>   |

**DIC**

Condenser filters:

- 10x/0,3 ----> **I**
- 20x/0,5 ----> **II**
- 40x/0,3 ----> **II**
- >40x ----> **III**

**Plus filter ring position "8" DIC**

**Phase contrast**

- Ph1 ----> 10x and 20x
- Ph2 ----> 40x
- Ph3 ----> 100x